



Product Code . DS-A-10029

Horizontal Autoclave

Description

Horizontal Autoclave

Features:

1. Sterilizing programs for fabric, instrument, rubber and liquid.
2. The positive pressure air-exhaust system to exhaust the cold air from chamber.
3. Build-in efficient steam generator, microcomputer temperature control.
4. Touch screen display of temperature, pressure, time, running status and failure warning.
5. Equipped with printer to print date, time and process parameters.
6. Post-drying function makes the residual humidity of sterilized articles less than 3%.
7. Perfect four-level authority management system to avoid mis-operation.
8. Over-heat auto protection device, dual door safety lock device, over-pressure relief valve, short-circuit protection and over-pressure protection.

Specifications:

Capacity	100L	135L	185L	300L

Chamber Size	400*850mm	500*710mm	500*950mm	632*1000mm
Chamber material	304 medical stainless steel			
Max working pressure	0.23Mpa			
Working temp	105~136℃			
Temp accuracy	±0.1℃			
Power consumption	11kw			17kw
Power supply	380V, 50HZ			
External size(mm)	1250*730*1660	1110*750*1742	1350*750*1742	1400*890*1742
Packaging size(mm)	1450*1000*1950	1370*890*2060	1510*890*2060	1700*1050*2060
Gross weight	310 kg	470kg	540 kg	830 kg

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/