



Product Code . DS-A-10031

Hand Wheel Vertical Pulse Vacuum Autoclave

Description

Hand Wheel Vertical Pulse Vacuum Autoclave

Features:

1. Pulse vacuum exhaust to ensure the effect of steam sterilization.
2. LCD screen display shows pressure, temperature, time and error codes.
3. Steam water internal circulation.
4. Removable probe detects the internal temperature of liquid directly, to ensure the sterilization effect of liquid program.
5. Standard built-in printer for recording the information during the entire sterilization process.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/