


Product Code . DS-A-10034

Table Top Autoclave Class S



Description

Table Top Autoclave Class S

Description:

As a type of high pressure sterilizer, it takes steam as its sterilization medium which is fast, safe and economic. They are common used in stomatology department, ophthalmology department, operating room and make the sterilization for wrapped or unwrapped instrument, fabric, utensils, culture medium, unsealed liquid, etc.



Features:

Build-in open type water tank

- a. Easy-clean open type water tank
- b. Water quality monitor
- c. Water filter inside
- d. Water level sensor
- e. A fully injected water tank can support repeated program running.

Specification:

- Volume: 24 L, 45 L, 80 L
- Design Pressure: -0.1/0.3MPa
- Design Temp.: 144?
- Temp. Range: 105~138?
- Chamber Material: Medical Stainless steel 304
- Standard Loading Tray: Type A, Type B
- Power Supply: 110/220V±10%, 50/60Hz, 220V±10%, 50/60Hz
- Power Consumption: 2.6 Kw, 5.8 Kw, 7.4 Kw
- Chamber Size(mm): ?250*450, ?316*618, ?358*717
- External Size(mm): 660*530*450, 830*640*550, 1002*730*560
- Packing Size(mm): 880*720*700, 1015*750*730, 1180*830*795
- Gross Weight(kg): 80, 140, 205.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/