

Product Code . DS-A-10035

Table Top Autoclave Class B



Description

Table Top Autoclave Class B

Features:

Build-in open type water tank

a. Easy-clean open type water tank

b. Water quality monitor

c. Water filter inside

d. Water level sensor

e. A fully injected water tank can support repeated program running.

Specifications:

Capacity	18 L	24 L	45 L	60 L	80 L
Design Pressure	-0.1/0.3MPa				
Design Temp.	144?				
Temp. Range	105~138?				
Chamber	Medical Stainless steel 304				

Material					
Standard Loading Tray	Type A		Type B		
Power Supply	110/220V±10%, 50/60Hz		220V±10%, 50/60Hz		
Power Consumption	2.6 Kw	2.6 Kw	5.8 Kw	7.4 Kw	7.4Kw
Chamber Size(mm)	250*350	250*450	316*618	358*595	358*717
External Size(mm)	560*530*450	660*530*450	830*640*550	1002*730*560	1002*730*560
Packing Size(mm)	780*720*680	880*720*700	1015*750*730	1180*830*795	1180*830*795
Gross Weight(kg)	70	80	140	185	205

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/