



Product Code . DS-A-10041

Table Top Dental Sterilizer

Description

Table Top Dental Sterilizer

Features:

- The newly designed type with opening water tank on the top, which makes washing and cleaning much more convenient.
- Microcomputer control, LCD display and simple operation interface.
- Small capacity design more suitable for clinic requirement, saving space and moving easily.

Technical Parameters:

Capacity	2L	8L
Sterilization Class	Class N	Class B
Sterilization Temp.	121?, 134?	
Chamber Size(mm)	188×175	170×320
Chamber Material	304 medical stainless steel	
Max Working Pressure	0.25Mpa	
Max Working Temp.	145?	
Power Supply	AC 110/220V, 50/60Hz	
Power	1.2KW	1.4KW

Consumption		
External Size(mm)	350*280*360	540*370*360
Packing Size(mm)	415*345*445	615*440*425
Gross Weight(kg)	36	45

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/