



Product Code . DS-A-10042

Hot Air Sterilizer

Description

Hot Air Sterilizer

Application:

Hot air sterilizer is used to sterilize small surgical instrument, glass, petri dishes etc.

This device is suitable for medical and veterinary clinics, hospitals, pharmacies, health care centers, and laboratories.

Features:

- Over-temperature alarm.
- Top air outlet with adjustable air volume.
- Adjustable sterilization time.
- Heating up fast, forced-air convection, fast sterilization.
- High precision microprocessor temperature control; 304 stainless steel inner chamber.
- Fan inside, fully sealed inside, no air circulation with external air, temperature rise fast.

Specifications:

Capacity	25L	50L	70L	105L	
Temp. Range	RT+10~250?				
Temp. Precision	0.1?				
Temp. Fluctuation	±0.5?				
Ambient Temp.	5~40?				
Timing Range	1~9999min				
Shelves No.	2 pcs				
Power Consumption	600W	900W	1150W	1600W	
Power Supply	220V/110V±10%,50/60Hz				
Internal Size (W*D*H) mm	300*300*270	420*350*350	400*375*500	500*350*500	500*350*500
External Size (W*D*H) mm	600*470*420	710*520*520	600*550*860	840*550*860	840*550*860
Packing Size (W*D*H) mm	690*550*600	840*680*630	680*680*990	970*680*990	970*680*990
Gross Weight(kg)	45	50	75	80	

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/