



Product Code . DS-A-10043

Portable Autoclave

Description

Portable Autoclave

Range of application:

It's desk-top autoclave sterilizer which is designed for the hospital, clinic, lab, etc. It's mostly suitable for sterilizing of surgical, dental instruments, glassware, culture medium and biological dressing, food and goods, etc.

Features:

- Easy operation, safe and reliable.
- Two safety valves keep double safety.
- It has self-adjustment function.
- It is well sealed and effectively prevents steam leakage.
- Electric heating or gas heating can be used.
- The sealing ring is made of silicon rubber.
- Over temperature & over pressure auto-protection.

-
- The sterilizer is provided with a power switch, user can control the operation easily.
 - Steam-water inner circulation system? No steam discharge and the environment for sterilizing will be clean and dry.
 - Immersion type electric heating pipe, thermal efficiency is high, stainless steel material tube is durable to use.
 - When at lower water level, or at the end of heat preservation time, have buzzer alarm, timing control type will automatically stop heating.
 - Double groove pressure gauge, could clearly indicate, convenient observation, reliable performance.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/