

**Product Code . DS-A-10044**

## Cassette Sterilizer



### Description

#### Cassette Sterilizer

##### Description:

It is used for the repaid sterilization of small instruments in stomatology, ophthalmology, E.N.T. department, gynecology, operating room, etc.

It only takes 6-7 minutes for the shortest sterilization cycle which can extend service life of the instrument.

Three shortcut keys for sterilization programs, the sterilization temperature and time can be set freely.

Small dimension, light weight, easy to operate, the cassette can be also used for instrument transfer and temporary storage.



Alarm for unqualified water and water shortage, the autoclave will stop working when there is a low temperature and low pressure.

##### Specifications:

Capacity	1.8L	5.2L	6.0L
Chamber Size(W*D*H)mm	280*180*38	380*180*78	480*180*78

External Size(W*D*H)mm	570*415*170	580*460*190
Temp. Range	115?~135?	
Temp. Fluctuation	±5?	
Ambient Temp.	5~40?	
Working Pressure	42 kPa ~212kPa	
Chamber Material	304 Stainless Steel	
Water Tank Volume	3.4 L	
Standard Accessories	Instrument loading tray*1, Waste water tank*1 Drain hose*1, Printing paper*2	
Power Supply	AC 220V/110V±10%, 50Hz/60Hz	
Packing Size(W*D*H)mm	580*580*700	
Net Weight(kg)	34	46
Gross Weight(kg)	49	61

---

Didac Scientific Pvt. Ltd,  
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
 Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)