



Product Code . DS-BAS-10048

Analytical Balances

Description

Analytical Balances



Description:

- Automatic internal adjustment in the case of a change in temperature $> 0,5\text{ }^{\circ}\text{C}$ and time-controlled every 4 hours
- High-stability mode and other filter settings can be selected
- Simple recipe weighing and documenting with a combined tare/print function. In addition, the ingredients for the recipe are numbered automatically and printed out with their corresponding number and nominal weight
- Identification number: 4 digits, printed on calibration protocol freely programmable
- Automatic data output to the PC/printer each time the balance is steady

Features:

- Fully automatic manufactured weighing cell from one piece of material
- Shock proof construction

-
- High corner load performance
 - Stable temperature behaviour
 - Short stabilisation time: Steady weight values within approx. 5 sec under laboratory conditions
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/