

Product Code . DS-BAS-10052

## Compact Portable Balance 200g



### Description

---

#### Compact Portable Balance 200g

##### Description:

Laboratory Equipment is the best choice for simple operation and economy. The tough durable ABS housing is designed to withstand field and laboratory environments, while being easy to clean and protected from accidental spills. With built-in ShockProtect protection Core balances can handle excessive overloads without a problem. Complete with a removable draft shield, AC adapter, plus integral security slot, Core balances are ready to work right out of the box.

##### Features:



- 9 weighing units
- Full range tare
- Zero Tracking
- Simple 4 button operation
- Backlit LCD display
- Stainless steel pan

- 
- Lightweight and portable
  - Low cost and simple to use
  - External calibration
  - Battery power (6 x AA Batteries)
  - Auto sleep / power down function to save battery life
  - AC Adapter

### **Specifications:**

- Capacity: 200g
- Readability: 0.01g
- Repeatability (S.D.): 0.01g
- Linearity ( $\pm$ ): 0.02g
- Pan Size: 4.7" Ø / 120mm Ø
- Display: Backlit LCD with 0.6" / 18mm high digits
- Material: ABS plastic
- Power Supply: 12 VAC, 150mA adapter or 6 x AA batteries
- Overall Dimensions: 6.8x10x3.4" / 173x255x86mm

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)