

Product Code . DS-BAS-10053

Compact Balance 200g



Description

Compact Balance 200g

Description:

Features for quick and easy weighing. Every Dune is equipped with an AC adapter, auto calibration through the keypad, low battery indicator, and overload indicator. With its alternate battery operation, splash-proof keypad, and sturdy plastic construction throughout, the Dune delivers accurate readings in the lab or in the field.

Features:

- 3 Weighing units (g, lb, oz)
- Full range tare
- Zero Tracking
- Large easy to read LCD display
- Overload protection
- Color coded, sealed keypad

-
- Lightweight and portable
 - Simple user-friendly operation
 - Precision load cell technology
 - Splash-proof to protect from accidental spills
 - External calibration
 - Battery power (6 x AA Batteries)
 - Low battery indication
 - Auto sleep / power down function to save battery life
 - AC Adapter

Specifications:

- Capacity: 200g
- Readability: 0.1g
- Repeatability (S.D.): 0.1g
- Linearity (\pm): 0.2g
- Pan Size: 5.7x5.7" / 145x145mm
- Units of Measurement: Gram (g), Ounce (oz)
- Calibration: External mass.
- Display: LCD with 0.6" / 16mm high digits
- Material: ABS plastic
- Power Supply: 12 VAC, 150mA adapter or 6 x AA batteries
- Overall Dimensions: 5.8x8.2x1.6" / 147x208x40mm