



**Product Code . DS-BAS-10054**

## Wash Down Scale 2000g

### Description

---

#### Wash Down Scale 2000g

##### Description:

Wash Down Scales are ideal for food preparation, food processing plants, agricultural facilities, and more, scales feature an wash down housing and a large stainless steel weighing pan. Color indicators for Under, Ok, and Over limits make check weighing fast and easy, and a brilliant backlit LCD with capacity tracker, indicators for overloads, low battery, and auto power-off, make the WBW very user friendly. They operate with an AC adaptor or rechargeable battery pack for portability; both are included.

##### Features:

- 5 weighing units (kg, g, lb, oz, lbs:oz)
- Full range tare
- Zero Tracking
- Check-weighing with low and high limits
- Check weighing LED bar graph

- 
- Large backlit LCD display
  - Capacity tracker
  - Large stainless steel pan
  - Non-slip adjustable levelling feet
  - Color coded, sealed keypad
  - Simple user-friendly operation
  - Precision load cell technology
  - Lightweight and portable
  - External calibration
  - Rechargeable battery
  - Auto sleep / power down function to save battery life
  - Low battery indication
  - AC Adapter



### **Specifications:**

- Capacity: 5 lb / 2000g
  - Readability: 0.0005lb / 0.2g
  - Repeatability (S.D.): 0.2g
  - Linearity ( $\pm$ ): 0.4g
  - Pan Size: 8.3x6.8" / 210x173mm
  - Units of Measurement: Gram (g), Kilogram (kg), Pound (lb), Pound-Ounce (lb-oz)
  - Calibration: External mass.
  - Display: Backlit LCD with capacity tracker
  - Material: ABS Plastic housing with Stainless Steel pan
  - Power Supply: 12 VDC, 800 mA adapter and internal rechargeable battery
-

---

- Overall Dimensions: 9.1x10.4x6" / 231x265x153 mm

---

Didac Scientific Pvt. Ltd,  
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)