

Product Code . DS-BAS-10056

Triple Beam Balance 610g



Description

Triple Beam Balance 610g

Description:

The Triple Beam Balance offers the functions and durability you need, plus an array of special features that you will appreciate. Every Triple Beam Balance comes with a below-balance weighing, an integral security bracket, a mass holder.

Features:

- Low cost and simple to use
- Magnetic damping for faster stability
- Ruler for measuring cm
- Weight holders for storing additional weights
- Rigid construction for durability
- Zero adjust facility to enable zero reading
- Stainless steel weighing pan

-
- Tripod legs for density
 - Large stainless steel pan
 - Lock down mounting slot for Kensington™ type lock and cable

Specifications:

- Class: Triple Beam
- Capacity: 610g
- Readability: 0.1g
- Repeatability (S.D.): 0.2g
- Linearity (±): 0.2g
- Pan Size: 6" Ø / 152mm Ø
- Units of Measurement: Gram (g)
- Calibration: External mass.
- Material: Aluminium
- Overall Dimensions: 19.9×4.3×6.3 inches / 505×110×160mm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/