

Product Code . DS-BAS-10059

Bench Counting Scale 4000kg



Description

Bench Counting Scale 4000kg

Description:

Laboratory Equipments full-featured CBD has the best value in an advanced parts counting scale. For the more complex counting applications, the flexible three displays and the ability to connect to a remote base. A secondary scale base, the CBD will meet all your parts counting needs.

Features:



- Internal counting resolution 1:400,000
- Preset tare function
- Full range tare
- Zero Tracking
- Memory accumulation of count
- Easy access to the rechargeable battery
- Stainless steel pan

-
- Overload protection
 - Pre-set counting with alarm
 - Programmable for counting to a pre-set number of parts
 - Large backlit LCD display
- b Full color coded keypad with numeric entry
- Precision load cell technology
 - Simple user-friendly operation
 - RS-232 bi-directional interface
 - External calibration
 - Rechargeable battery
 - Auto sleep / power down function to save battery life
 - Low battery indication
 - AC Adapter
 - Supports remote scale of up to 4 load-cells for enhanced weighing and counting
 - Non-slip adjustable leveling feet

Specifications:

- Capacity: 8 lb / 4000 kg
 - Readability: 0.0002 lb / 0.1 g
 - Repeatability (S.D.): 0.0002lb / 0.0001kg
 - Linearity (\pm): 0.0004lb / 2g
 - Pan Size: 8.9x10.8" / 225x275mm
 - Units of Measurement: Pound (lb), Kilogram (kg)
 - Interface: Bi-directional RS-232
 - Calibration: External mass.
-

-
- Display: LCD, backlit
 - Material: ABS plastic, stainless steel platform
 - Power Supply: Main power cord / Internal rechargeable battery
 - Overall Dimensions: 12.4×14×4.3 inches / 315×355×110mm
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/