



Product Code . DS-CW-10062

Sets of Calibration Weights

Description

Sets of Calibration Weights

Description:

Complete sets of stainless steel calibration weights of standards, with weights ranging from 1mg to 200g. Each set is housed in a rugged lined aluminium storage case, and supplied complete with glove and tweezers.

Size: Set of 23 Stainless Steel Calibration Weights in Case 1mg – 200g, Set of 11 Stainless Steel Calibration Weights in Case 1g – 200g.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/