



Product Code . DS-C-10066

Mini Centrifuge Machine

Description

Mini Centrifuge Machine

Features:

- Timing function
- Rotor with 8 positions
- Convex type switch, easy to use
- Long life, reliable motor
- Completely transparent cover
- Low noise, stable running
- Humanization design, switch control for cover
- Unique rotor clip design, easy to install and replace
- Large LED display, indicate speed and set time
- Auto stop working when cover opens during running and restart running after cover closed.

Parameters:

Max Speed	4000rpm	6000rpm	7000rpm	10000rpm	Adjust Max 1
Max RCF	1350xg	3000xg	6000xg	7500xg	Adjust
Rotor	8*1.5ml/0.5ml/0.2ml				
Noise	?45dB		?50dB		
Power Supply	AC110/220V±10%?50/60Hz				
External Size	175*156*125				
(W*D*H) mm					
Package Size	260*200*150				
(W*D*H) mm					
Net/Gross Weight(kg)	1.4/3				

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
🌐 www.didacscientific.co.uk/