

**Product Code . DS-C-10067**

## Large Capacity High Speed Refrigerated Centrifuge



### Description

#### Large Capacity High Speed Refrigerated Centrifuge

#### Features:



- \* LCD display
- \* Brushless motor
- \* Over 100 selectable procedures and freely self-defining modes.
- \* Safety: over-speed and over-temperature protection, imbalance detection, shock absorber to keep motor running smoothly, three layers steel structure to secure personnel and machine safety.

#### Specifications:

Max. Capacity	4*250ml
Max. Speed	20600rpm
Max. RCF	29400xg
Standard Rotor	12*1.5/2.0ml
Speed Precision	±50rpm
Refrigeration System	Import CFC-free compressor and control valve (R404a)
Temp. Control Range	-20°C~40°C
Temp. Precision	±1°C
Control and Driver System	Frequency motor, micro control, direct driver

Timing Range	0?99h59min
Noise	<65dB
Identify Rotor	Automatic
Consumption	1.8KW
Power Supply	AC220V, 50Hz
External Size(W*D*H) mm	630*760*400
Packing Size(W*D*H) mm	890*890*585
Net/Gross Weight(kg)	140/162

---

Didac Scientific Pvt. Ltd,  
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
 Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)