



Product Code . DS-C-10069

Table Top Refrigerated Centrifuge

Description

Table Top Refrigerated Centrifuge

Features:

- LCD display: indicating setting parameters and actual parameters.
- Stainless steel chamber.
- Inching function.
- Mute mechatronics motor door lock.
- Brushless motor: high torque and variable frequency.
- Safety: over-speed and over-temperature protection, imbalance detection, shock absorber to keep motor running smoothly.
- Over 20 selectable procedures, 9 levels speed-up and speed-down control.

Specifications:

Max. Capacity	4*100ml		
Max. Speed	16000rpm	20000rpm	
Max. RCF	18180xg	27700xg	

Standard Rotor	12*1.5/2.0ml
Speed Precision	±50rpm
Refrigeration System	Import CFC-free compressor and control valve (R404a)
Temp. Range	-20°C?40°C
Temp. Precision	±1°C
Timing Range	0?99h59min
Noise	<65dB
Consumption	1.5KW
Power Supply	AC220V, 50Hz
External Size(W*D*H) mm	710*630*350
Packing Size(W*D*H) mm	860*780*550
Net/Gross Weight(kg)	73/95

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/