

Product Code . DS-C-10072

Refrigerated Centrifuge



Description

Refrigerated Centrifuge

Features:

- Microprocessor controls , LED displays.
- Maintenance-free brush less motor.
- 20 standard programs, programmability is available.
- Automatic rotor recognition and imbalance detection.
- Adjustable acceleration levels: 15, deceleration levels: 15.
- Automatic temperature compensation.
- Double Lid and Lock system, door locks automatically once starting running.
- The shortest Acc/Dec time is less than 15s.
- Pre-cooling and fast cooling function, temperature dropped to 40? within 10 minutes.

