



Product Code . DS-C-10074

Table Top High Speed Centrifuge

Description

Table Top High Speed Centrifuge

Features:

1. DC Brush less motor with high start torque and frequency inversion for smooth running, maintenance-free.
2. External cold-rolled steel with powder coating steel, stainless steel centrifugal chamber.
3. Microprocessor control: LED display , Memory flashing function, conversion of RPM/RCF.
4. Procedures of adjustable speed from grad 0 to 9 and 40 kinds of self-defining models.
5. Functions: Electronic lock lid, imbalance detection, over-speed protection, inching function, short-spin function.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details ✉ sales@didacscientific.co.uk
💻 www.didacscientific.co.uk/