



Product Code . DS-C-10075

## Large Capacity Refrigerated Centrifuge

### Description

---

#### Large Capacity Refrigerated Centrifuge

##### Features:



- LCD display, microprocessor control, indicating setting parameters and actual parameters at the same time.
- Imported drive system and automatic control system to ensure steady and reliable operation.
- Safety: Over-temperature and over-speed protection, electronic locking, imbalance detection.
- AV motor: high torque and variable frequency, free from contamination caused by carbon dust.
- Can set 16 procedures, parameters can be modified during running.
- 10 levels speed-up and speed-down control.
- Can indicate rotor speed ascending and declining curves, centrifugal integral curves.

##### Specifications:

Max. Capacity	6*1000ml
Max. Speed	10000rpm
Max. RCF	18300xg
Speed Precision	±50rpm

Control and Driven System	Frequency motor, micro control, direct driver
Refrigeration System	Import CFC-free compressor and control valve (R404a)
Temp. Range	-20°C~40°C
Temp. Precision	±1°C
Timing Range	0?99h59min
Noise	<68dB
Consumption	5.7KW
Power Supply	AC220V, 50Hz
External Size(W*D*H) mm	840*730*950
Packing Size(W*D*H) mm	1085*835*1165
Net/Gross Weight(kg)	280/360

---

Didac Scientific Pvt. Ltd,  
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
 Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)