

Product Code . DS-C-10076

PCR Centrifuge



Description

PCR Centrifuge

Features:

- Microprocessor control?LED display.
- Easy to set and read RCF?speed synchronizes in time.
- Brushless motor, free maintenance, speed up and down quickly.
- Electric lid lock, over speed and imbalance protection.

Specifications:

Max. Capacity	5ml*10
Max. Speed	16000rpm
Max. RCF	19040*g
Speed Accuracy	±20rpm
Temp. Range	-20?~40?
Temp. Accuracy	±1.0?
Timing Range	1~99min
Noise	?55dB
Power Supply	AC110/220V±10%, 50/60Hz
External Size(L*W*H) mm	480*480*275

Packing Size(L*W*H) mm	630*560*450
Net Weight(kg)	43
Gross Weight(kg)	57

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/