



Product Code . DS-C-10077

## Crude Oil-water Test Centrifuge

### Description

---

#### Crude Oil-water Test Centrifuge

##### Features:



- Microprocessor control
- LED display: speed, RCF, time, and temperature, conversion of RPM/RCF.
- External cold-rolled steel with powder coating, stainless steel centrifugal chamber and lid.
- DC Brushless motor, special damping device make the motor run steady, free-maintenance.
- Alarm: Lid open, imbalance, over speed, rotor mounting abnormally.
- Function: automatically electronic lid lock, over speed, over temperature and imbalance protection, emergency lid-lock release, inching function, short-spin function, Memory flashing function.

##### Specifications:

|                |         |
|----------------|---------|
| Max. Capacity  | 4*100ml |
| Max. Speed     | 5000rpm |
| Max. RCF       | 4500xg  |
| Speed Accuracy | ±30rpm  |
| Temp. Range    | 10?80°C |
|                |         |

|                         |              |
|-------------------------|--------------|
| Temp. Accuracy          | ±2°C         |
| Timing Range            | 0?99min59s   |
| Noise                   | <58dB        |
| Consumption             | 1.75KW       |
| Power Supply            | AC220V, 50Hz |
| External Size(W*D*H) mm | 610*720*1140 |
| Packing Size(W*D*H) mm  | 760*870*1340 |
| Net Weight(kg)          | 108          |
| Gross Weight(kg)        | 123          |

---

Didac Scientific Pvt. Ltd,  
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
 Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)