



Product Code . DS-C-10078

Automatic Uncovering Centrifuge

Description

Automatic Uncovering Centrifuge

Features:

- LCD display, microprocessor control.
- Stainless steel chamber.
- Can store 20 procedures.
- 10 levels speed-up and speed-down control.
- Mute mechatronics motor door lock.
- Brushless motor: variable frequency, free from contamination caused by carbon dust.

Specifications:

Capacity	48/76*7/5ml	124*7/5ml	124*7/5ml
Max. Speed	4000rpm	5000rpm	5500rpm
Max. RCF	378?xg	378?xg	480?xg
Speed Precision	±50rpm		
Refrigeration System	/		
			Import CFC-free compressor and control

			valve (R404a)
Temp. Range	/		-20oC?40oC
Temp. Precision	/		±1°C
Control and Drive System	Frequency motor, micro control, direct driver		
Timing Range	0?99h59min		
Noise	<65dB		
Consumption	800W	1000W	2000W
Power Supply	AC220V, 50Hz		
External Size(W*D*H) mm	586*644*370	710*560*770	710*560*770
Packing Size(W*D*H) mm	780*840*570	900*750*970	900*750*970
Net/Gross Weight(kg)	50/75	75/105	150/200

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/