

Product Code . DS-C-10079

Table Top Dairy Centrifuge



Description

Table Top Dairy Centrifuge

Features:

1. Microprocessor control.
2. The digital display indicates the speed, time and temperature.
3. Stainless steel centrifugal chamber and lid.
4. Electric lid lock, over speed, over temperature and imbalance protection.

Specifications:

Max. Capacity	8*Gerber tubes	4*300ml	4*800ml
Max. Speed	1500rpm	5000rpm	
Max. RCF	504xg	4730xg	
Speed Accuracy	±50rpm	±20rpm	
Temp. Range	RT+10?~65?		RT~90?
Timing Range	0~99h59min	0?99min	
Noise	?60dB	?58dB	
Consumption	1.5KW	600W	450W
Power Supply	AC220V, 50Hz	AC110V/220V±10%,50/60Hz	
External Size(L*W*H)	520*415*330mm	620*450*360mm	685*500*395mm

Packing Size(L*W*H)	720*610*530mm	820*520*610mm	870*580*610mm
Net Weight	35kg	52kg	67kg
Gross Weight	60kg	68kg	77.6kg

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/