



Product Code . DS-C-10081

Hematocrit Centrifuge

Description

Hematocrit Centrifuge

Features:

- * Microprocessor control
- * Stainless steel chamber.
- * Maintenance-free brushless DC motor.

Specifications:

Max Speed	12000rpm
Max RCF	13800xg
Max Capacity	12*10ml
Timer Range	0~99h59min
Noise	<65dB
Power Supply	AC110/220V, 50/60Hz
External Size(W*D*H)mm	500*440*330
Package Size(W*D*H)mm	750*640*530
Net Weight(kg)	30
Gross Weight(kg)	50

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/