



Product Code . DS-C-10082

Blood Bag Centrifuge

Description

Blood Bag Centrifuge

Features:

- * Microprocessor control, digital display with the temperature, time and speed.
- * Stainless steel centrifugal chamber and lid.
- * Brushless motor, free maintenance, speed up and down quickly.
- * Automatically electric lid lock, over speed, over temperature and imbalance protection.

Specifications:

Max. Capacity	1000ml*4	1200ml*6	2400ml*6
Max. Speed	6000rpm	6000rpm	7000rpm
Max. RCF	6600xg	6600xg	11900xg
Speed Accuracy	±20rpm		
Temp. Range	-20°C ~ 40°C		
Temp. Accuracy	±1°C		
Refrigeration Function	Yes		
Timing Range	1~99h59min		
Noise	<55dB	<65dB	<60dB
Power Supply	AC110V/220V±10%, 50/60Hz		380V, 50/60Hz

External Size(L*W*H) mm	800*675*830	750*850*940	940*890*1000
Packing Size(L*W*H) mm	920*800*1320	1000*880*1350	1200*1050*1565
Net/Gross Weight(kg)	260/305	260/310	620/685

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/