



Product Code . DS-C-10084

PRP Centrifuge

Description



PRP Centrifuge

Features:

- * Microprocessor control, LED display: rotor model, time, RCF and speed.
- * Double Lid and Lock system.
- * RPM setting function, and short-spin function.
- * Automatic rotor recognition and imbalance detection?audible alarm for error and imbalance.
- * Easily adjustable speed, RCF, time and temperature during the running without stop.
- * Acc/Dec: choice of 15 acceleration and deceleration rates.
- * Quiet, maintenance-free brushless drive.
- * 20 programs memory, programmability is available.
- * Easily setting and reading RCF, speed synchronize in time Double layers steel structure, spray anti-corrosion coating outside.
- * The shortest Acc/Dec time is less than 28s.
- * Built-in exhaust filtering devices, it can prevent aerosol pollution during the centrifugal process.

Specifications:

Max. Speed	5000rpm
Max. RCF	5000xg
Speed Accuracy	±10rpm
Noise	<62dB
Timing Range	1?9999min, inching function
External Size(W*D*H) mm	425*510*415
Net Weight(kg)	50
Gross Weight(kg)	67

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/