



Product Code . DS-C-10086

Micro High Speed Centrifuge

Description

Micro High Speed Centrifuge

Features:

- * LCD display: indicating setting parameters and actual parameters, including RCF.
- * Stainless steel chamber.
- * Variable frequency brushless motor.
- * Can store 20 procedures.
- * Mute mechatronics motor door lock.
- * 10 levels speed-up and speed-down control.

Specifications:

Max. Capacity	12*10ml
Max. Speed	16500rpm
Max. RCF	18780xg
Timing Range	0-99h59min
Control and Drive System	Frequency motor, micro control, direct drive
Noise	<65dB

Consumption	350W
Power Supply	AC220V, 50Hz
External Size(W*D*H) mm	380*260*310
Net/Gross Weight(kg)	20/28

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/