



**Product Code . DS-C-10094**

## **Mini Centrifuge Single Speed 4000rpm**

### **Description**

---

#### **Mini Centrifuge Single Speed 4000rpm**

##### **Description:**

Compact in design and personal sized, Mini Centrifuge Single Speed 4000rpm mini-4 centrifuge is ideal for micro filtration and quick spin downs from the walls and caps of micro centrifuge tubes.

This Mini Centrifuge Single Speed 4000rpm Mini-4 is supplied with rotors and adapters to accommodate 0.2ml 0.5ml 1.5ml 2.0ml tubes as well as 0.2 mL strips and individual tubes. equipped with 2 types of rotors and several tube carriers for 2.0ml / 1.5ml / 0.5ml / 0.2ml tubes and 0.2ml x 8 PCR tube rack. The rotors are easily interchanged.

1. LED display the rolling speed and timing.
2. The rotor is clip-on. It is very convenient for rotors exchanging and pick-and-place.
3. The cover switch is outward designed and superstructure automatically opens 95 degrees which is convenient for one-hand operation.
4. Silent running and stable. free-maintenance motor, long life and safe.
5. Large radius rotor with 8 tube holes, centrifugal force is 1.5 times bigger than ordinary 6-hole rotor at the same rolling speed.
6. Delicate appearance. Multi-purpose to meet different experiments demands.

---

## **Applications:**

- 1.Extract serum from blood
2. Extract the supernatant from a variety of samples
3. Rapid subsidence of samples
4. Separation of trace blood cells
5. Microbial sample processing
6. PCR experimental division centrifugal

## **Specifications:**

- Max Speed: 4000rpm
- Max RCF: 850g
- Max. Sample Capacity: 8x2.0ml/1.5ml and 8x0.5ml and 2x8x0.2ml
- Time Range: 1s ~9999min or continuous
- Safety: self-lock, over-speed, over-temp., status diagnosis system
- Power: 45W
- Speed Range: 6000rpm
- Acceleration to Max. Speed: 15 seconds
- Deceleration to Stop: 15 seconds
- Rotor: 8x2.0ml/1.5ml centrifuge tubes
- Strip Rotor: 2 PCR strips or 2x8x0.2ml tubes
- Tube adapter: 8x0.5ml and 8x0.2ml