



Product Code . DS-C-10095

Low Speed Centrifuge (Economical Type)

Description



Low Speed Centrifuge (Economical Type)

Features:

An apparatus consisting essentially of a compartment spun about a central axis to separate contained materials of different specific gravities, or to separate colloidal particles suspended in a liquid.

Specifications:

- Max Speed: 15,000
 - Safety: self-lock, over-speed, over-temp., status diagnosis system
 - Other Functions: speed/RCF switch, touch operation, operating status display, alert
 - Power: 105W
 - Speed Range: 100 ~15000rpm (increment: 100rpm)
 - Acceleration to Max. Speed: 15 seconds
 - Deceleration to Stop: 15 seconds
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/