

Product Code . DS-C-10097

Table Top Centrifuge



Description

Table Top Centrifuge

Features:

- * Brush less motor
- * Stainless steel chamber
- * Can store 20 operation procedures.

Specifications:

Max. Capacity	8*20ml	6*50ml	4*100ml	6*50ml	4*250ml	4*100ml
Max. Speed	4000rpm		5000rpm		6000rpm	
Max. RCF	2200xg		4420xg	4020xg	5140xg	
Timing Range	0~99h59min	1~60min	0~99h59min	0~99h	0~99h59min	
Speed Precision	±30rpm	/	±50rpm	±30rpm	±50rpm	
Noise	<60dB	<65dB		<58dB	<65dB	
Consumption	100W	180W	350W	300W	500W	
Power Supply	AC220V, 50/60Hz					
External	380*260*310	395*340*305	330*295*255	600*540*360		

Size(W*D*H) mm					
Net Weight(kg)	18	20	25	20	35
Gross Weight(kg)	33	35	44	28	50

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/