

Product Code . DS-C-10099

Table Top Low Speed Refrigerated Centrifuge



Description

Table Top Low Speed Refrigerated Centrifuge

Features:

- * LCD display, intelligent control, touchable and press button.
- * Stainless steel chamber.
- * Automatically record the parameters, can start by RCF directly.
- * Mute mechatronics motor door lock.
- * Brush less motor: high torque and variable frequency, combined with shock absorbers.
- * Over 20 selectable procedures and freely self-defining modes.
- * 10 levels speed-up and speed-down control.

Specifications:

Max. Capacity	4*300ml		
Max. Speed	5500rpm	6500rpm	8000rpm
Max. RCF	4800xg	4800xg	5940xg

Speed Precision	±50rpm
Refrigeration System	Import CFC-free compressor and control valve
Temp. Range	-20oC-40oC
Temp. Precision	±1°C
Control and Drive System	Frequency motor, micro control, direct drive
Timing Range	0~99h59min
Noise	<65dB
Consumption	1.8KW
Power Supply	AC220V, 50Hz
External Size(W*D*H) mm	630*760*400
Net/Gross Weight(kg)	110/135

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/