

Product Code . DS-C-10100

Low Speed Refrigerated Centrifuges



Description

Low Speed Refrigerated Centrifuges

Features:

- * Type: Table-Top
- * LCD display.
- * Stainless steel chamber.
- * Maintenance-free brushless DC motor.
- * Automatically electric lid lock, over speed, over temperature, and imbalance protection.

Specifications:

Max Speed	5000rpm
Max RCF	4390xg
Max Capacity	4×250ml
Speed Accuracy	±50rpm
Time Range	LCD: 1s?99h59min59s
Temp. Range	-20???40?
Temp. Accuracy	±1?
Refrigerant	CFC-free refrigerant

Noise	?65dB
Power Supply	AC110/220V, 50/60Hz
External Size (W×D×H) mm	640*680*380
Net Weight/ kg	95
Gross Weight/ kg	120

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/