

Product Code . DS-C-10101

Large Capacity High Speed Refrigerated Centrifuge DID202



Description

Large Capacity High Speed Refrigerated Centrifuge DID202

Features:

- * Microprocessor control
- * Brush less motor

Specifications:

Max. Capacity	4*750ml	4*250ml
Max. Speed	21000rpm	20600rpm
Max. RCF	30700xg	29400xg
Speed Accuracy	±30rpm	±50rpm
Temp. Range	-20~40°C	
Temp. Accuracy	±1°C	
Timing Range	0~99h59min59s	0~99h59min
Noise	<58dB	<65dB
Go? sumption	2.5KW	1.8KW
Power Supply	AC220V, 50Hz	
External Size(W*D*H) mm	800*700*400	630*760*400
Net Weight(kg)	110	140
Gross Weight(kg)	125	162

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/