



Product Code . DS-C-10102

High Speed Centrifuge Machine

Description

High Speed Centrifuge Machine

Features:

- * LCD display: indicating setting parameters and actual parameters, including RCF.
- * Stainless steel chamber.
- * Variable frequency brushless motor.
- * Can store 20 procedures.
- * Mute mechatronics motor door lock.

Specifications:

Max. Capacity	6*50ml	4*100ml	12 x 10ml	6*100ml
Max. Speed	16500rpm	18500rpm	16600rpm	18500rpm
Max. RCF	18780xg	23610xg	19200xg	26019xg
Timing Range	0~99h59min	0~99h59min59s	1s~999min	
Noise	<65dB		<60dB	
Consumption	350W	450W	600W	800W
Power Supply	AC220V, 50Hz			
External	500*440*330	400*350*340	330x390x325	380*450*390

Size(W*D*H) mm				
Net/Gross Weight(kg)	30/50	20/28	30/38	42/50

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/