



Product Code . DS-C-10103

Table Top Refrigerated Centrifuge Machine

Description

Table Top Refrigerated Centrifuge Machine

Features:

- * Brush less motor
- * 10 levels speed-up and speed-down control.
- * Over 20 selectable procedures and freely self-defining modes
- * Acc/Dec: choice of 15 acceleration and deceleration rates.

Specifications:

Max Capacity	300ml*4		
Max Speed	5500rpm	6500rpm	8000rpm
Max RCF	5400xg	4800xg	5940xg
Speed Precision	±10rpm		±50rpm
Refrigeration System	/ Import CFC-free compressor and control valve		
Temp. Range	-20?~40?		
Temp. Precision	±1?		
Control and drive system	/ Frequency motor, micro control, direct drive		
Timing Range	1s~999min/continuous/s 0~99h59min		

	hort spin	
Noise	?60dB	?65dB
Consumption	1.8kw	
Power Supply	AC220V, 50Hz	
External Size(W*D*H)mm	450*575*520	630*760*400
Net/Gross Weight(kg)	110/145	110/135

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/