



Product Code . DS-C-10104

Low Speed Refrigerated Centrifuge Machine

Description

Low Speed Refrigerated Centrifuge Machine

Features:

- * Brush less motor, microprocessor control.
- * CFC-free refrigeration, fast cooling, cycle defrost system.
- * External cold-rolled steel with powder coating steel, stainless steel centrifugal chamber and lid.
- * 10 procedures of adjustable speed from 0 to 9 and 40 kinds of self-defining modes.
- * 20 selectable procedures and freely self-defining modes, 10 levels speed-up and speed-down control.

Specifications:

Max. Capacity	4*750ml	4*750ml	4*750ml	6*1000ml
Max. Speed	5500rpm	6500rpm	8000rpm	6000rpm
Max. RCF	4800xg	4800xg	5940xg	6680xg
Speed Precision	±50rpm			±20rpm
Refrigeration System	CFC-free(R404a)			CFC-free (R404a or R134a)
Temp. Range	-20?~40?			

Temp. Precision	$\pm 1^{\circ}\text{C}$	
Timing Range	0~99h59min	0~59h59min59s
Noise	<65dB	<58dB
Consumption	2.0kw	4.0kw
Power Supply	AC220V, 50/60Hz	
External Size(W*D*H) mm	710*560*770	830*720*1250
Net/Gross Weight(kg)	160/208	268*318

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/