



Product Code . DS-C-10105

Table Top Low Speed Large Capacity Centrifuge

Description

Table Top Low Speed Large Capacity Centrifuge



Features:

- * LCD display: speed, RCF, time, conversion of RPM/RCF
- * External cold-rolled steel with powder coating, stainless steel centrifugal chamber and lid
- * Programmable microprocessor control (AC brushless motor).
- * Electronic lock lid and interlock safety system. Auto lock rotor system
- * Fan-forced air circulation system for fast recovery of temperature
- * 10 procedures of adjustable speed from levels 0 to 9 and 40 kinds of self-defining modes

Specifications:

| | |
|----------------|---------------|
| Max. Capacity | 4*750ml |
| Max. Speed | 5500rpm |
| Max. RCF | 5310xg |
| Speed Accuracy | ±20rpm |
| Timing Range | 0~99h59min59s |
| Noise | <58dB |
| | |

| | |
|-------------------------|--------------|
| Consumption | 2.5KW |
| Power Supply | AC220V, 50Hz |
| External Size(W*D*H) mm | 800*700*400 |
| Net Weight(kg) | 110 |
| Gross Weight(kg) | 135 |

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/