

Product Code . DS-FD-10112

Vertical Freeze Dryer





Description

Vertical Freeze Dryer

Description:

Vertical Freeze Dryer widely used in drugs, biological products, chemical and food industries. On the heat-sensitive substances such as antibiotics, vaccines, blood products, hormones and other biological tissue enzymes, freeze-drying technology is applied.

1. LCD touch screen.
 2. With historical data inquiry function.
 3. With USB interface to export the data.
 4. Drying chamber is transparent which is visual and safe.
 5. With cascade refrigeration technology.
 6. Ice condenser trap and operation panel are made by stainless steel.
 7. Large capacity ice condenser trap and without coil inside.
 8. Pre-freezing shelf can be used as guide barrel to speed up the drying rate.
 9. Shelves are stainless steel, which can be adjustable per the requirements.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/