



Product Code . DS-FHASC-10116



Class II A2 Biological Safety Cabinet

Description

Class II A2 Biological Safety Cabinet

Advantages:

1. Motorized front window.
 2. Large LCD display, all information displayed.
 3. Automatic air speeds adjustable with filter block.
 4. Single piece SS removable work tray.
 5. Side & back wall is made of single piece stainless steel
 6. Front 100 slanted to offer operator comfort while working for long time, reduce glare and maximize reach in to work area.
 7. Interlock function: UV lamp and front window; UV lamp and blower, fluorescent lamp; blower and front window.
 8. With memory function in case of power-failure.
 9. Armrest is above the intake of grill to enable farther reach inside the cabinet without hampering safe airflows inside the cabinet working area.
 10. Differential pressure sensor for monitoring of condition of HEPA filters.
-

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/