



Product Code . DS-FHASC-10119

Electric Height Adjustable Base Stand DID206

Description

Electric Height Adjustable Base Stand DID206

Features:

1. Heavy load and smart balance
2. Low maintenance cost
3. Concise and easy with nice appearance
4. Easy to install and disassemble
5. Smooth running speed and low noise
6. Max Adjustable Route 400 mm
7. Multipurpose for laboratory apparatus
8. Customize as you want

Specifications:

Height of Base Stand	Min. Height(mm)	480	580	680	730	780

	Middle height(mm)	530	680	730	905	980
	Max. Height(mm)	580	780	980	1080	1180
Max Load		600kg (100-600kg optional)				
Load speed		3mm/s				
Working frequency		10%, max continuous working time is 2 minutes				
Noise		<42dB				
Work Temperature		+5?~40?				
Material		Pedestal: Square steel pipe(40*40*2mm)				
		Lifting Rod :Aluminium alloy				
Power Supply		AC220V±10%,50/60Hz;110V±10% 60Hz				
Controller Consumption		100~200W				
Caster		Foot master Caster				

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details ✉ sales@didacscientific.co.uk
 🌐 www.didacscientific.co.uk/