



Product Code . DS-FHASC-10121

Walk-in Fume Hood

Description

Walk-in Fume Hood



Advantages:

1. Operators are available to walk into the fume hood to operate if need.
2. Large room exhaust system, suitable for large laboratory instruments and equipment.

Specifications:

External Size (W*D*H)	1200*800*2200mm	1500*800*2200mm	1800*800*2200mm
Internal Size (W*D*H)	980*600*1650mm	1280*600*1650mm	1580*600*1650mm
Max Opening	1300mm		
Exhaust Airflow Volume	1150m ³ /h	1500m ³ /h	1850m ³ /h
Inflow Air Velocity	0.3~0.8m/s		
Noise	62dB		
Vibration	XYZ direct ?4?m		
Illumination	?800Lux		
Fluorescent Lamp	12W*2	16W*2	16W*2
Power Supply	AC220V±10%, 50/60Hz;110 V±10%, 60Hz		
Consumption	400W	600W	800W
Front Window	Manual, Toughened glass		
Material	Exterior	1.0mm cold-rolled steel with bacteria powder coating.	
	Interior	High grade melamine board with good acid and alkali	

		resistance function
--	--	---------------------

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/