



Product Code . DS-FHASC-10124

Ducted Fume Hood

Description

Ducted Fume Hood

Advantages:

1. Adjustable air speed: 9 levels.
2. Manual front glass window, height adjustable.
3. With memory function in case of power-failure.
4. Resistant to moderate acid and alkali.
5. PP Built-in axial flow blower, melamine board inward, resistant to acid, alkali, and anti-corrosion.

Specifications:

External Size(W*D*H)	1000*800*2365 mm	1200*800*2365 mm	1500*800*2465 mm	1800*800*2465 mm
Internal Size(W*D*H)	790*635*820mm	990*635*820mm	1290*635*820 mm	1590*635*820 mm
Work Surface Height	900mm			
Max Opening	750mm			
Air Velocity	0.3~0.8m/s			

Noise	?60dB			
Exhaust Duct	PVC, Standard length:4 meters			
	?300mm			
LED Lamp	12W *1	12W *1	16W *1	12W *2
Blower	Built-in PP axial flow blower; Speed adjustable		Built-in PP centrifugal blower; Speed adjustable	
Front Window	Manual, 5mm toughened glass, height adjustable			
Power Supply	AC220V±10%,50/60Hz;110V±10%, 60Hz			
Consumption	400 W	400 W	500 W	500 W

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/