



Product Code . DS-FHASC-10128

Ductless Fume Hood

Description

Ductless Fume Hood

Advantages:

1. LCD Touch screen control panel, easy to operate.
2. With memory function in case of power-failure
3. Double-layer structure: 1mm sheet metal surface, Chemical resistant phenolic resin work table.
4. Temperature and humidity sensors, can detect indoor temperature and humidity.
5. Electronic control system, anti-overload, anti-electric shock, stable performance, long service life.
6. 8°slope frontergonomics design, fatigue-free working posture.
7. Three side transparent Acrylic windows, front window reversal design, easy to operate.
8. Inside and outside probe, detect indoor air pollution and filter conditions. Audible and visual alarm for changing filter.

Specifications:

External Size	1000*800*2120mm	1200*800*2120mm	1500*800*2120mm	1800*800*2120mm
---------------	-----------------	-----------------	-----------------	-----------------

(W*D*H)				
Internal Size (W*D*H)	910*690*705mm	1110*690*705mm	1410*690*705mm	1710*690*705mm
Work Surface Height	900mm			
Max. Opening	650mm			
Air Velocity	0.3~0.6m/s			
Airflow Volume	220m3/h		440m3/h	
Noise	?55dB(A)			
Fluorescent Lamp	11W*2	11W*2	11W*3	11W*4
UV Lamp	Emission of 253.7 nano meters for most efficient decontamination			
Blower	Built- in centrifugal blower, speed adjustable			
Front Window	Acrylic window; Manual; reversal design;			
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz			
Consumption	400W	400W	500W	500W

Didac Scientific Pvt. Ltd,
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
 Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/