



**Product Code . DS-FHASC-10129**

## Ducted Fume Hood DID210

### Description

#### Ducted Fume Hood DID210

#### Advantages:



1. Microprocessor control system, LED display.
2. Built-in centrifugal blower.
3. With memory function in case of power-failure.
4. Folding acrylic front window, down part with free stop function.
5. Interlock function: UV lamp and blower, fluorescent lamp.

#### Specifications:

External Size(W*D*H)	700*620*1150mm
Internal Size(W*D*H)	640*550*700mm
Max Opening	660mm
Air Velocity	0.3~0.8m/s
Noise	?60dB?A?
HEPA Filter	99.999% efficiency at 0.3um
LED Fluorescent Lamp	4W *2
Blower	Built-in centrifugal blower; Speed adjustable

Front Window	Folding acrylic window, down part with free stop function
Consumption	200W
Power Supply	AC220V±10%,50/60 Hz; 110V±10%,60 Hz

---

Didac Scientific Pvt. Ltd,  
 Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom  
 Direct Contact Details  [sales@didacscientific.co.uk](mailto:sales@didacscientific.co.uk)  
 [www.didacscientific.co.uk/](http://www.didacscientific.co.uk/)