



Product Code . DS-FHASC-10131

A2 Biological Safety Cabinet

Description

A2 Biological Safety Cabinet

Description:

Three protection: operator, sample and environment. Airflow system: 70% air recirculation, 30% air exhaust A2 Cabinet is suitable for work with microbiological research in the absence of volatile or toxic chemicals and radionuclide.



Advantages:

1. Motorized Front Window.
2. Time Reserve Function.
3. Automatic air speed adjustable with filter block.
4. Work area surrounded by negative pressure, it can ensure maximum safety in work area.
5. HEPA Filter life and UV Life indicator.
6. Most accessories are standard. No need to pay more, save your money.
7. Audio and visual alarm (Filter replacement, Window over height, abnormal Air flow velocity etc).

8. Interlock Function: UV lamp and front window; UV lamp and blower, fluorescent lamp; Blower and front window.

9. Foot Switch. Adjust front window height by foot during experiment, to avoid airflow turbulence caused by arm movement.

10. Remote Control. All functions can be realized with it, making the operation much easier and more convenient.

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/