

Product Code . DS-FHASC-10132

Class III Biological Safety Cabinet



Description

Class III Biological Safety Cabinet

Description:

Class III Biosafety Cabinet Totally enclosed and gas-tight with ULPA filtered supply and exhaust air. Work is performed with long-sleeved gloves. The cabinet is kept under negative pressure of at least 120 pa, and airflow is maintained by a dedicated exterior exhaust system. It can protect the operator, product and environment.

Advantages:

1. UV lamp ensures work area sterilization effectively.
2. An angled cabinet front ensures an ergonomic working posture.
3. Exhaust air is double-filtered through high-quality ULPA filters with typical efficiency of 99.999% for 0.1~0.2µm particles