



Product Code . DS-FHASC-10134

Necessary Tools for Filter Replacement

Description

Necessary Tools for Filter Replacement

1. Formalin Fumigation Sterilizer

- Sterilize the cabinet with formalin for 8 hours
- Power: 180w
- Volume: 100ml
- Power supply: 220V/50Hz

2. Ammonium Hydrogen Carbonate Neutralizer

- Power supply: 220v/50Hz
- Power: 180W
- Volume: 100ml
- Neutralize the balance formalin with ammonium hydrogen carbonate

3. Protective Garment

4. Protective Gloves

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/