



Product Code . DS-FHASC-10135

Biosafety Accessories

Description

Biosafety Accessories

Infrared Sterilizer

Infrared sterilizer is a gasless, flameless sterilizer. It can prevent dangerous spattering of microorganisms by the infrared heat during the operations of sterilizing the contaminated loops, needles, mouth of tube and small bottle.

Specifications:

The maximum temperature is towards 930° in 10 minutes after turned on.

Two types:

- Type I, 37mm dimension of cavity
- Type II, 20mm dimension of cavity
- Size: 150x95x210mm
- Weight: 1600g
- Central high-temperature: 825°±50°
- Voltage: 220 V/50Hz



-
- Power: 180W

Armrest

- Alleviate fatigue on operator's arms
- Adjustable height range: 85~100mm
- Prevent air obstruction

Laboratory Chairs

- Economy line- artificial leather
- Washable, easy to clean, robust and durable
- Anatomically designed seat and backrest
- Backrest adjustable in height and inclination
- Seat height adjustable by gas lift
- Height range: 38-48cm

Didac Scientific Pvt. Ltd,
Martinfield Business Centre, 108 Martinfield, Welwyn Garden City , United Kingdom
Direct Contact Details  sales@didacscientific.co.uk
 www.didacscientific.co.uk/